## Human Th1/Th2 Staining Kit

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This package insert must be read entirely before using this product. For proper performance, use the insert provided with each individual product received.

**Catalog Number**

KTH101 - 25

KTH101 - 100

Optimization for human Th 1/Th2 staining in anticoagulated blood.

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**For research use only. Not for use in diagnostic procedures.**

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## INTRODUCTION

The kit is designed to detect frequencies of Th 1 and Th2 cells in T lymphocytes (or CD4<sup>+</sup> helper T cells) in human anticoagulated blood. Th 1, Th2 and Th 17 cells are subsets differentiated from helper T cells (strictly speaking, that is Th0 cells) under physiological and pathological conditions. There are extremely few Th 1, Th2 or Th 17 cells in peripheral blood in a resting state (unstimulated state, e.g. normal physiological state in human) due to the weak differentiation ability of Th0 cells, thus, it's hardly to detect IFN- $\gamma$ , IL-4 and IL-17A. Th0 cells could differentiate to Th 1, Th2 or Th 17 cells dependent on different cytokines when stimulated by external factors (e.g. stimuli, pathogens), and more IFN- $\gamma$ , IL-4 and IL-17A could be detected. In experiment, it is to measure the ability of Th cells in response to stimuli.

Generally, Phorbol 12-Myristate 13-Acetate (PMA) and Ionomycin are used as stimuli. PMA, a PKC activator, could activate phosphorylation of many protein kinases in downstream, cascade reaction induces protein expression, then results in activation of T cells. PKC could be activated in the combination of DAG and Ca<sup>2+</sup> in cells, with the involvement of Ionomycin (a calcium transporter, transport Ca<sup>2+</sup> from organelles to cytoplasm), PKC could be further activated in T cells. So PMA and Ionomycin could activate T cells synergistically.

Activated T cells would secrete cytokines to extracellular environment, which could not be detected by flowcytometry, thus cytokines should be blocked in the cell. Cytokines are synthesized in the Golgi and some may traffic through the endoplasmic reticulum to be released in soluble form in the extracellular environment. Blocking secretion by disrupting Golgi to cut off the way of cytokine transport. Generally, Brefeldin A and/or Monensin are used as blocker.

Besides, how to choose Th cells by correct gating is important. The surface marker of Th cells is CD3<sup>+</sup>CD4<sup>+</sup>, however, CD4<sup>+</sup> cells will decrease obviously, even disappear after PMA stimulation in human anticoagulant blood for 4 hours due to endocytosis of CD4 molecules induced by PMA. Generally, we use CD3 and CD8 to gate CD4<sup>+</sup> cells, that is, CD3<sup>+</sup>CD8<sup>-</sup> cells would be regarded as CD3<sup>+</sup>CD4<sup>+</sup> cells.

The choice of anticoagulant would influence the assay. It will not affect activation of Th cells and cytokine secretion when using heparin sodium or heparin lithium as anticoagulant in anticoagulated blood, while using EDTA or sodium citrate as anticoagulant, no cytokines could be detected. In this case, PBMCs should be used for stimulation rather than anticoagulated blood, however, the frequency of cytokines detected may be lower than that in heparin-anticoagulated blood.

## LIMITATIONS OF THE PRODUCT



The product is intended for flowcytometry applications.

The product should not be used beyond the expiration date on the label.


Do not mix reagents with those from other lots or sources.

## Related products/相关产品

Cat No.	Product Name	Size
CS0001	Phorbol 12-Myristate 13-Acetate (PMA)	100 $\mu$ l
CS0002	Ionomycin Calcium	50 $\mu$ l
CS0003	Brefeldin A (BFA)	60 $\mu$ l
CS0004	Monensin	100 $\mu$ l
F0001	多聚甲醛溶液, 4%	100 ml
GAS005/2	FIX&PERM Kit	100T
KTH001	Human Th1/Th2/Th17 Staining Kit	25T/ 100T
KTH117	Human Th 17 staining kit	25T/ 100T
KTH201	Mouse Th 1/Th2 staining kit	25T/ 100T
KTH217	Mouse Th 17 staining kit	25T/ 100T
KTR101	Human regulatory T Cell staining kit	25T/ 100T
KTR201	Mouse regulatory T Cell staining kit	25T/ 100T
LSB01	Lysing solution for FACS 10 $\times$	100T
LSC01	FCM Lysing solution for BC (ready-to-use)	100T
LYS01	FCM Lysing solution (Fixative Free) 10 $\times$	100T

## Examples of results/结果示例

**Flow cytometric analysis of Human Th1/Th2 Staining Kit.** The staining pattern of IFN- $\gamma$  and IL-4 on resting heparin anticoagulated blood (A<sub>0</sub>, B<sub>0</sub>) and PBMCs from EDTA anticoagulated blood (C<sub>0</sub>, D<sub>0</sub>). PMA/Ionomycin stimulated heparin anticoagulated blood (A<sub>1</sub>, B<sub>1</sub>) and PBMCs from EDTA anticoagulated blood (C<sub>1</sub>, D<sub>1</sub>). Dot plot analysis are derived from gated CD3<sup>+</sup>/CD8<sup>+</sup> cells population. Flow cytometry was performed on a Thermo Fisher Attune NxT. The health status of the volunteers is unclear.

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
**Note:** PMA/Ionomycin Mixture (250×) and BFA/Monensin Mixture (250×) can be stored at - 20°C, FIX & PERM Medium A and Medium B can be stored at room temperature, other reagents can be stored at 2 - 8°C. All reagents can be stable for at least 1 year when stored at recommended condition.

**OTHER SUPPLIES REQUIRED**

1. **Medium** (e.g. RPMI 1640, DMEM).
2. **Fetal Bovine Serum** (Maybe needed)
3. **Compensation Beads** (Maybe needed)
4. **Lymphocyte separation medium** (Maybe needed)
5. **Paraformaldehyde, 4%** (Cat No: F0001, MultiSciences) (Maybe needed)
6. **12 × 75 mm round bottom test tubes.**
7. **Vortexer.**
8. **CO<sub>2</sub> incubator.**
9. **Swing-out horizon centrifuge** (with rotor for 15 ml tubes).

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- Add 1 ml 1× Flow Cytometry Staining Buffer to each tube. Centrifuge at 300 × g for 5 minutes, discard supernant.
- Resuspend pellet by adding 500 μl 1× Flow Cytometry Staining Buffer to each tube, ready for determination. Or add 500 μl 1 - 4 % paraformaldehyde to resuspend pellet, then store at 2 - 8℃ , protect from light, determine within 24 hours.

#### Detection by flowcytometry

- Voltage and compensation adjustment

It is recommended to adjust voltage and compensation of flow cytometer by using compensation beads. It is important to adjust voltage of FSC/SSC and gate target cells because compensation beads are often smaller than cells.


- Correct gating to determine the frequencies of Th 1/Th2 cells in CD3<sup>+</sup>CD8<sup>+</sup> T cells.

**Note:** Acquire at least 20,000 to 30,000 CD3<sup>+</sup>CD8<sup>+</sup> T cells. Depending on the donor, frequencies of cytokine producing cells derived from activation of human PBMCs can vary widely for a particular cytokine. In order to make statistically significant frequency measurements, sufficiently large sample sizes should be acquired during flow cytometric analysis.

In particular, the number of IL-4 producing cells can be very low or even negligible on PMA/Ionomycin stimulated cells. In these cases, IL-4 polarization cultures should be considered.

#### Excitation wavelength and emission wavelength of fluorophore

Fluorophore	Ex (nm)	Em (nm)
FITC	495	519
R-Phycoerythrin (PE)	480;565	578
PerCP-Cy5.5	488	682
APC	650	660

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## PROTOCOL

### Sample collection

Collect appropriate amount of anticoagulant blood by using heparin as anticoagulant or other anticoagulants (e.g. EDTA, sodium citrate), store at room temperature or 2 - 8°C , and detect within 24 hours.

### Sample preparation

- 1a. *For anticoagulant blood with heparin as anticoagulant*, pipet 125  $\mu$ l anticoagulant blood into 12  $\times$  75 mm round bottom test tubes, add 125  $\mu$ l medium without FBS and 1  $\mu$ l PMA/Ionomycin Mixture (250 $\times$ ) and 1  $\mu$ l BFA/Monensin Mixture (250 $\times$ ). Mix 125  $\mu$ l anticoagulant blood and 125  $\mu$ l medium, served as control. Mix well, incubate in CO<sub>2</sub> incubator at 37°C for 4 - 6 hours, vortex every 1 - 2 hours during incubation.
- 1b. *For anticoagulant blood with other anticoagulant (e.g. EDTA, sodium citrate)*, isolate peripheral blood mononuclear cells (PBMCs) by using Lymphocyte separation medium . Resuspend pellet at 1  $\times$  10<sup>7</sup>/ml in medium with 10 % FBS. Pipet 250  $\mu$ l PBMCs into 12  $\times$  75 mm round bottom test tubes, add 1  $\mu$ l PMA/Ionomycin Mixture (250 $\times$ ) and 1  $\mu$ l BFA/Monensin Mixture (250 $\times$ ). With PBMCs only as control. Mix well, incubate in CO<sub>2</sub> incubator at 37°C for 4 - 6 hours, vortex every 1 - 2 hours during incubation.

**Note:** PMA/Ionomycin Mixture (250 $\times$ ) and BFA/Monensin Mixture (250 $\times$ ) are extremely volatile, tighten immediately the lid after use.

2. Pipet 100  $\mu$ l sample or control into new test tubes, add 5  $\mu$ l Anti-Human CD3, FITC and 5  $\mu$ l Anti-Human CD8 $\alpha$ , PerCP-Cy5.5 to each tube. Vortex to mix well, and incubate at room temperature for 15 minutes, protect from light.
3. Add 100  $\mu$ l FIX & PERM Medium A to each tube. Vortex to mix well, and incubate at room temperature for 15 minutes, protect from light.
4. Add 1 ml pre-cool 1 $\times$  Flow Cytometry Staining Buffer to each tube. Centrifuge at 300  $\times$  g for 5 minutes, discard supernant.

**Note:** Discard supernant as much as possible, with little residual liquid left.

5. Add 100  $\mu$ l FIX & PERM Medium B, 5  $\mu$ l Anti-Human IFN- $\gamma$ , PE and 5  $\mu$ l Anti-Human IL-4, APC to each test tube. Vortex to mix well, and incubate at room temperature for 15 minutes, protect from light.

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